



B10 DATA EXPORT PROCESS FLOW

The documentation for the export format consists of information on the structure and syntax for how the EHI will be exported by the product such as, for example, Consolidated-Clinical Document Architecture (C-CDA) document(s) or data dictionary for comma separated values (csv) file(s).

Navigate to Administration>>Practice Setup>> Data Export.



To Export Patient Data, Click on Data Export>>Click on +New Export.



Once Clicked on +New Export button the below screen will be displayed.

The screenshot shows the 'Data Export' setup screen. At the top, there is a navigation bar with 'Dashboard' and 'Practice Setup' tabs. Below this is a menu bar with options: General, Billing, Schedule, EMR, Allergy, Meaningful Use, Menu Privileges, Password Expiration, Register IP Address, and Data Export (highlighted with a red box). The 'Data Export' sub-menu is also highlighted with a red box. The main content area has a 'Select All Patients' checkbox. Below it are fields for 'Patient Name' (a dropdown menu), 'Archive Name' (a text input), 'File Type' (a dropdown menu), and 'Encryption Key' (a text input). There is a 'Comments' section with a dropdown menu showing 'Tahoma' and a rich text editor with various formatting icons. At the bottom, there are 'Submit' and 'Cancel' buttons.

The Requestor has the option to select the Patient by clicking on the Patient Name option, Has the ability to select the Single or Multiple Patient by Name or Account No# as shown below.

This screenshot shows the same 'Data Export' setup screen, but with a patient selection modal open. The 'Patient Name' dropdown now shows the number '2'. The modal has an 'OK' button and a '1 Selected' indicator. It contains a table with two columns: 'Account #' and 'Name'. The table has one row with the values '2' and 'Priyanka-july 's, Priyanka-june'. There are checkboxes in the first column of the table. At the bottom of the modal, there is a pagination bar showing 'Page 1 of 1' and 'Record(s) 1 - 1 of 1'. The 'Submit' and 'Cancel' buttons are still visible at the bottom of the main screen.

Account #	Name
2	Priyanka-july 's, Priyanka-june

Multiple Patient Selection Screen as shown below.

The screenshot shows the 'Data Export' window with the 'Multiple' patient selection method. A table lists one patient record with account number 4 and name 'namtsae, evitcanI'. The 'Select All Patients' checkbox is unchecked.

Account #	Name
4	namtsae, evitcanI

Requestors also have the option to select all the patients by just clicking the check box of Select All Patients as shown below.

The screenshot shows the 'Data Export' window with the 'Select All Patients' checkbox checked. The patient selection table is no longer visible, and the 'Archive Name' and 'File Type' fields are empty.

Other Options and functionalities:



Archive Name: To Name the File

File Type:

Encryption Key: To secure the Data files with key/Password.

Comments: Comments box to add any comments or notes or message.

Submit Button: To Submit the request.



Clicking on the submit button, will create the request for the selected Patients or Multiple Patients or All Patients.

A backend program will be scheduled every day to export the requested data. The status can be seen on the Data Export Screen.

This program will be performed/completed in 3 steps:

Step 1: The request is raised/submitted, and the status will be as New on Data Export Screen.

Step 2: The requested data is once generated in folders; the status will be changed to InProgress.

Step 3: Once the data is copied from physical server to Cloud environment, a new URL will be generated, and the Status will be as Completed.

By using this URL, the requestor can download the files.

The URL(Downloads) has the expiry date which can be seen on the screen.

Note: The Download option is available only for the Requestor (who raised the request), other users don't see the download option.

On Data Export Screen we see the Archive Name, Date Generated, Created By, Status and Expiry Date.

Step 1: The request is raised, and the status will be as New on Data Export Screen.

Data Export				
New Export				
Archive Name	Date Generated	Created By	Status	Expiry Date
DevAsthma7_130	11/20/2023 2:11:08 AM	Nagamalla Suman	New	
DevAsthma7_129	11/20/2023 2:10:42 AM	Nagamalla Suman	New	

Step 2: The requested data is generated in folders, the status as InProgress.

Data Export				
New Export				
Archive Name	Date Generated	Created By	Status	Expiry Date
DevAsthma7_127	11/15/2023 5:31:18 AM	Nagamalla Suman	Inprogress	
DevAsthma7_126	11/14/2023 7:35:33 AM	Nagamalla Suman	Inprogress	
DevAsthma7_125	11/14/2023 7:35:13 AM	Nagamalla Suman	Inprogress	
DevAsthma7_124	11/14/2023 7:34:36 AM	Nagamalla Suman	Inprogress	
DevAsthma7_123	11/3/2023 9:43:49 AM	admins modulemd	Inprogress	
DevAsthma7_122	11/3/2023 7:25:59 AM	admins modulemd	Inprogress	


Step 3: Status as Completed and with Expiry Date and Download option.

The Download option/icon is available/visible only for the Requestor as shown below.

Data Export				
New Export				
Archive Name	Date Generated	Created By	Status	Expiry Date
DevAsthma7_130	11/20/2023 2:11:08 AM	Nagamalla Suman	Completed	12/20/2023 12:00:00 AM
DevAsthma7_129	11/20/2023 2:10:42 AM	Nagamalla Suman	Completed	12/20/2023 12:00:00 AM
DevAsthma7_128	11/19/2023 7:40:40 AM	Indra Chauhan	Completed	12/19/2023 12:00:00 AM


Download

Guest/other User view as shown below, as mentioned above the other User cannot see the Download option/icon.




Archive Name	Date Generated	Created By	Status	Expiry Date
DevAsthma7_130	11/20/2023 2:11:08 AM	Nagamalla Suman	Completed	12/20/2023 12:00:00 AM
DevAsthma7_129	11/20/2023 2:10:42 AM	Nagamalla Suman	Completed	12/20/2023 12:00:00 AM
DevAsthma7_128	11/19/2023 7:40:40 AM	Indra Chauhan	Completed	12/19/2023 12:00:00 AM

Click on Download option to download the Data Files as shown below:




Download the Folder and Select the Patient folder to view C-CDA in XML files as shown below.



Imp Note: Single Patients can download the C-CDA Files from the Patient Portal.

Login to Patient Portal>> Documents.



Under Documents>> File Name, C_CDA Hyperlink to view/download the data.

