



## Export Format & Naming Conventions

### A. General Format Notes

- CDA for both single-patient and patient population **Electronic Health Information (EHI)** is exported in **C-CDA R2.1** format.
- Exported data is determined by selected conditions.
- If no checkboxes are selected, only Patient Data is exported.
- If checkboxes are selected, exported files vary accordingly

#### 1. Default Export Location:

📁 [\\SERVERNAME\BarcodeScans\HL7ExportFiles\CDAexports\\_#](#) (# = Session ID)

| Name               | Date modified     | Type                 | Size |
|--------------------|-------------------|----------------------|------|
| CDAXMLExport_9.Zip | 2/19/2025 2:43 PM | Compressed (zipp...) |      |

#### 2. XMLDATAExport\_#.zip

📁 [\\SERVERNAME\BarcodeScans\HL7ExportFiles\CDAexports\\_#\XMLDATAExport\\_#.Zip](#)

##### Contains

- **CHART\_DOCS;DataTable.xml** – contains all chart docs details for the exported date and Medical Record Number (MR) range
- **CHART\_STORE;DataTable.xml** – contains all chart store details for the exported date and MR number range
- **PATIENT;DataTable.xml** – contains all patient data table details for the exported date and MR number range
  - All patient data tables included (except those hidden via **HIDE** option)
  - See [Including and Excluding Data](#)

| Name                        | Type     | Compressed size | Password ... | Size  |
|-----------------------------|----------|-----------------|--------------|-------|
| CHART_DOCS;DataTable.xml    | XML File | 2 KB            | No           | 10 KB |
| CHART_STORE;DataTable.xml   | XML File | 1 KB            | No           | 7 KB  |
| PATIENT;DataTable.xml       | XML File | 4 KB            | No           | 39 KB |
| PT_ALLERGIES;DataTable.xml  | XML File | 2 KB            | No           | 11 KB |
| PT_COMPLAINTS;DataTable.xml | XML File | 1 KB            | No           | 6 KB  |

#### 3. CDAXMLExport\_1.zip

📁 [\\SERVERNAME\BarcodeScans\HL7ExportFiles\CDAexports\\_578\XMLDATAExport\\_#.Zip](#)

- Contains CDA files
- Included only with **Bulk CCDA Export**

| Temp1_CDAXMLExport_9.Zip ▶ 303735,9,9,129,CDASUMDATA.Zip |          |                 |          |
|----------------------------------------------------------|----------|-----------------|----------|
| Name                                                     | Type     | Compressed size | Password |
| 303735,9,9,129,CDAFULSUM.xml                             | XML File | 6 KB            | No       |

## B. Naming Formats

### 1. Chart Documents (CDAXMLDOC.Zip)

\\SERVERNAME\BarcodeScans\HL7ExportFiles\CDAexports\_578\CDAXMLExport\_#.Zip

- Contains chart document XML files
  - Documents exported as embedded PDFs in CCDA format
    - The fields in GREEN come from the Patient table, pa\_clinic and pa\_account fields
    - The fields in Pink are from chart\_store table cs\_mr\_num and cs\_cl\_key
    - The field in Red is the CD\_DOC\_ID field from the cdview\_cddocid field of the pt\_chart\_docs data table
- 247806,9,9,73,PT\_CHART\_DOCS;DataTable.xml
- Format:** [MRN, PCM LOGON CLINIC, PATIENT-CLINIC, PATIENT-ACCOUNT, Chart DOC ID, CDADOCID.xml]
  - (Example: 12345,1,1,2,55555,CDADOCID.xml)

| BarcodeScans ▶ HL7ExportFiles ▶ CDAexports_775 ▶ 286715,9,9,88,CDAXMLDOC.Zip |          |                 |              |        | Search 28.. |
|------------------------------------------------------------------------------|----------|-----------------|--------------|--------|-------------|
| Name                                                                         | Type     | Compressed size | Password ... | Size   |             |
| 286715,9,9,88,204505,CDADOCID.xml                                            | XML File | 7 KB            | No           | 16 KB  |             |
| 286715,9,9,88,204540,CDADOCID.xml                                            | XML File | 132 KB          | No           | 185 KB |             |
| 286715,9,9,88,206308,CDADOCID.xml                                            | XML File | 76 KB           | No           | 107 KB |             |

### 2. Patient Data Tables (CDAXMLDOC.Zip)

\\SERVERNAME\BarcodeScans\HL7ExportFiles\CDAexports\_1234\#,#,#,#,CDAXMLDOC.Zip

Multiple XML files per patient, one for each Patient Data Table.

Format: [MRN,CLINIC,ACCT-CLINIC,ACCOUNT-NUM,Data Table Name;DataTable.xml]

(Example: 12345,1,1,1,PT\_ALLERGIES;DataTable.xml)

## C. Data Export for Migrating to a Different EHR

When transitioning to a new Electronic Health Record (EHR) system, the new system must be able to import the exported data. The data tables are saved in Microsoft Data XML format, which allows them to be reloaded into a recordset for processing.

### 1. Understanding the XML Data Files

The exported data is stored as XML files using Microsoft ADO (ActiveX Data Objects) in Visual Basic 6.0. These files maintain the integrity of the database tables and can be reloaded for further processing in the new EHR system.

## 2. Example of How the Data is Saved

Below is a Visual Basic 6 example of how a data table is saved as an XML recordset document:

```
Dim FilePath as String
Dim TableName as String
Dim rstTemp As ADODB.Recordset
```

```
' Save the recordset as an XML file
rs.Save FilePath & "\" & TableName & ";DataTable.xml", adPersistXML
```

This process ensures that the data is stored in a structured XML format that can be reloaded into a recordset later.

## 3. Reloading the Data for Processing

Once saved, the data can be reloaded into a recordset for processing in the new EHR system using the following method:

```
Set rstTemp = New ADODB.Recordset
rstTemp.CursorLocation = adUseClient
rstTemp.CursorType = adOpenStatic
rstTemp.LockType = adLockBatchOptimistic
rstTemp.Open FilePath & "\" & TableName & ";DataTable.xml"
```

```
' Output the number of records in the recordset
debug.print rstTemp.RecordCount
```

This approach ensures that the data can be accessed and manipulated within the new EHR environment efficiently.

## 4. Next Steps

To successfully migrate data, ensure that:

- The new EHR system supports importing XML recordsets.
- The file structure and data integrity are maintained during export and import.
- Appropriate testing is conducted to validate the imported data before going live.

If you experience any issues during the export process, please contact [support@primeclinical.com](mailto:support@primeclinical.com)